


Less-invasive nasal test



The Coronavirus Ag Rapid Test Cassette (Swab) from MediPort Diagnostics

is an in vitro immunochromatographic assay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct nasal swab specimens directly from individuals who are both symptomatic and asymptomatic.

It is intended to aid in the rapid diagnosis of SARS-CoV-2 infections. Negative results from patients with symptom onset beyond ten days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.

The Coronavirus Ag Rapid Test Cassette (Swab) does not differentiate between SARS-CoV and SARS-CoV-2.

The kit contains 20 individual test kits, ready to use. Contents of the kit are:

- 20 Test Cassettes
- 20 Extraction Tubes with Buffer and Tips
- 20 Sterile Swabs
- 1 Workstation
- 1 Package Insert

Specification

Storage Condition: 2°C–30°C
 Sample type: nasal
 Shelf Life: 18 months

A CE approved product

Benefits of the MediPort COVID-19 Antigen Rapid Test Kit

- **UK Government tested & approved**
 Also proven to detect both VOC1 (Kent, UK) and VOC2 (South Africa) COVID-19 variants.
- **Quick**
 Get a positive or negative result within 15 minutes.
- **Reliable**
 97.25% sensitivity, 100% specificity and an overall accuracy of 98.73%.
- **Simple method**
 Less-invasive nasal 'Self Swab' sample collection.
- **Simple process**
 Testing that can be carried out or supervised by a trained person (Non-Medical).
- **Flexible**
 Tests both asymptomatic and symptomatic individuals
- **Affordable**
 A cost-effective solution for rapid testing
- **Scalable**
 We can supply in large volumes for efficient testing.

IVD Connect™ App



IVD Connect™ is a purpose-built end-to-end digital platform for all your organisations diagnostics and screening. Offering secure reader connectivity, it transforms traditional Lateral Flow IVD tests into web connected diagnostics with accessible, shareable data. It provides businesses a unique digital data management tool featuring:

End user data collection app

Collect test data & improve the experience for your end users who are testing at work or at home.

Live data & results tracking

Monitor data trends, results and outbreaks in real-time to enhance response time & disease control.

Scalable data management

Capture, store and manage your reader-generated data for research & commercial purposes

New features and functionality at no extra cost

Free updates and no ongoing service charges

Future integration with 3rd Party technologies

Vaccine Certificates, digital passports and NHS track and trace

Provided free for companies with 250+ employees

Provides all the features and benefits our IVD Connect digital platform with no upfront or ongoing costs to your organisation.

Pay as you go testing

No subscription or weekly fees at just £0.25 per test for companies less than 250 employees.

Benefits of the IVD Connect™ App

Providing a digital platform for fast, easy and affordable self-COVID-19 testing on both Android and IOS mobile platforms

Free for all business with 250+ employees when purchased with our products

Cloud Data Capture

via a recorded Photo of Test Result

Client Dashboard for managing and monitoring testing

- Easy to navigate user-friendly interface
- User Group Authentication
- Instantly share data with relevant stakeholders

Map Visualisation

- Geotagged, time-stamped results and interactive maps
- Monitor test results in different locations

Data Confidence & Security

All data collected is fully Secure, GDPR compliant and is data fully controlled by your organisation.

New features Coming soon

- Multi Site Dashboard Functionality
- Digital Test Result Interpretation (DTRI)
- Auto PCR in App Booking & Testing
- Vaccine and Digital Passport integration
- Fit to Fly Certificates

Train, Test & Track programme for businesses



Train

MediPort Diagnostics offers a unique COVID-19 test training solution for your workforce. We train your key staff to become COVID-19 test experts in delivering workplace testing, who then can supervise and train you employees to be confident in conducting self testing at home.

Our training can be delivered online or face to face by our UKAS accredited trainers.

Test

With fully competent and trained test specialists within your company you can now conduct supervised work place testing and employee self test training without the need costly medical or healthcare professional.

We also offer a dedicated COVID-19 Telephone support line just in case your staff need extra support from our dedicated medical team.

Track

Our **IVD Connect™** digital platform provides your staff simple test user experience to record their test via their mobile phone. Test data is then automatically transferred to on your own secure cloud based account. Managing your test all test data is simple with your own company dashboard and the export tool allows you to share your all your test results with PHE.

Key Benefits

- Makes staff testing more affordable and scalable
- Training model is more cost effective than managed testing services.
- Keep your workspace safe by training your staff to become more confident in testing themselves.
- MediPort diagnostic tests are easy to use with simple self swab sampling with no medical professionals needed to administer tests
- IVD Connect provides a complete digital management for your staff and organisation testing needs.
- Confidence that all data collected is fully secure, GDPR compliant and is data fully controlled by your organisation.
- To find out more about our train, test and track programme please get in touch.

Email: info@mediport.com

Visit: www.mediport.com